

PFL Switch Mount

Unique capabilities:

1. The PFL switch mount is very adjustable, rigid and durable. It can place switches in an unlimited number of positions and hold them securely without moving.
2. The PFL base clamp can be attached to a wheelchair easily without having to customize the base clamp. It can be attached to round tubing, non-round tubing, directly through the existing holes on the wheelchair, and hard mounted to the seat pan.
3. The PFL switch mount can mount switches, switch plates, environmental controls, cell phones, and other small light weight equipment.
4. The PFL switch mount can attach to wheelchairs, walkers, beds and desks.
5. The PFL switch mount can be mounted in a vertical or horizontal position and various other angles in between without impacting stability.
6. The PFL switch mount can be mounted on either the left or right side without having to purchase any adapters.
7. The PFL switch mount can be attached to the Profiler. Allowing easy positioning of the device and switch.
8. The PFL switch mount can swing away with the single lift of a lever.
9. The PFL switch mount is easy to remove either by opening the quick release lever or removing the pull pin to release all of the tubing from the base.
10. The PFL switch mount comes with all the tools and directions necessary to attach and position the mount.

Measurements & Options:

11. The base rotates 360 degrees and pivots out 270 degrees.
12. Standard adjustable lengths: 9 to 16 inches and 16 to 21 inches.
13. The compound angle switch head allows you to achieve numerous angles.
14. The footprint of the base clamp is only 2 ¼ x 2 ¼ inches. There are base clamp adapters for 1 inch square tubing and ¾, 7/8, and 1 inch round tubing.

Durability facts:

15. Aircraft Aluminum is used in the construction of the PFL switch mount making it durable, light weight and easy for teachers, family and care providers to work with.