

CJT Adjustable Table Stand

Key Facts:

- ❑ Great for assessment as it can be adjusted to accommodate the specific needs of numerous people.
- ❑ The CJT Table Stand can be assembled and adjusted without tools.
- ❑ The front of a device can be rotated 360 degrees of the device's center and tilted at a 45 degree angle up to a 45 degree angle down.
- ❑ With the unique flip rod adaptation a device can be positioned from just above the surface of a table to 17 inches above the surface of a table while maintaining stability.
- ❑ The steel, flat base design keeps the device stable at any adjusted position.
- ❑ The flat base legs allow paper work and books to be placed directly under the communication device.
- ❑ The CJT Table Stand can support all augmentative communication devices, especially eye gaze controlled devices.
- ❑ Perfect to use in conjunction with the Profiler, Profiler Angular, and ProLite mounting systems.

Overview:

The CJT Adjustable Table Stand is the optimum choice when it comes to table stands. Its innovative design allows you to achieve a variety of angles, heights and rotated positions without using tools.